

## FACT SHEET

# TRlactive® Select heel

Supporting and cushioning orthotic blank  
with heel spur relief

**Product features:**

- Thin TRlactive Technology in the cushion layer
- Yellow marking of the functional zone in cases of chronic heel and ankle pain
- Requires little modification and saves time
- Functional zones support, guide, and cushion the foot without noticeable transitions

**TRlactive Select heel**

**Item number:** 311 5005 0000 0\*\*

**Size\*\*** 35–48

**Treatment options:****Recommendations for use:**

- Chronic heel and ankle pain
- Heel spur
- Plantar fasciitis
- Arch decreased / splay foot
- Metatarsalgia

With chronic heel and ankle pain, TRlactive Select heel helps after a medical retailer has fitted it to the patient. The orthotic blank features a core with a cushioning cut-out in the heel area and an additional fan design to relieve the plantar fascia. In the cushion layer, materials featuring different densities are connected without causing irritation, thus supporting physiological gait.



More information is available  
in the online shop:  
[experts.bauerfeind.com](https://experts.bauerfeind.com)