

## FACT SHEET

# TRlactive® Select rigidus

Supporting and cushioning orthotic blank  
with rigidus spring

### Product features:

- Thin TRlactive Technology in the cushion layer
- Supporting core with cushioning base layer
- Rigidus spring made of thermoplastic fiberglass; available in the following versions:  
**reinforced on the left, reinforced on the right,**  
and reinforced on **both sides**
- Yellow marking of the indication-specific functional zone in cases of problems with the metatarsophalangeal joint
- Functional zones support, guide, and cushion the foot without noticeable transitions
- Requires little modification and saves time
- Reinforced heel edge stabilizes the foot

### TRlactive Select rigidus

**Item number** reinforced on the left: 312 6005 0000 0\*\*

**Item number** reinforced on the right: 312 7005 0000 0\*\*

**Item number** reinforced on both sides: 312 5005 0000 0\*\*


**Size\*\*** 35–48

### Treatment options:



### Recommendations for use:

- Hallux rigidus
- Hallux limitus
- Planovalgus foot
- Arch decreased / splay foot
- Metatarsalgia



In cases of hallux rigidus / hallux limitus, the TRlactive Select rigidus foot orthosis provides help after a medical retailer has fitted it to the patient. The orthotic blank combines a supporting, fiberglass-reinforced rigidus spring core and functional zone cushion layer. This orthotic blank is particularly suitable for acute hallux rigidus.



More information is available  
in the online shop:  
[experts.bauerfeind.com](https://experts.bauerfeind.com)