

## FACT SHEET

# TRlactive® Select

### Product features:

- **Medium-wide** and **wide** version
- Thin TRlactive Technology in the cushion layer
- Functional zones support, guide, and cushion the foot without noticeable transitions
- Requires little modification and saves time
- Reinforced heel edge stabilizes the foot

### TRlactive Select

**Item number** wide version: 313 8505 0000 0\*\*

**Item number** medium-wide version: 313 5005 0000 0\*\*

**Size\*\*** 35–48

### Treatment options:



### Recommendations for use:

- Planovalgus foot
- Arch decreased / splay foot
- Metatarsalgia
- Collapsed high-arched foot



Flat orthotic blank, requiring little modification, with a supporting core, sturdy EVA base layer and 3 mm absorbing functional zone cushion layer. Thanks to its thin design, the TRlactive Select is suitable for many shoe types.



More information is available  
in the online shop:  
[experts.bauerfeind.com](https://experts.bauerfeind.com)